

EASYLIFE

INSTRUCTIONS FOR CLEANING AND DISINFECTION

Caution

- For Clinical Use Only.
- The following processes are validated for the EASYLIFE nasal mask when used in accordance with the manufacturer of the cleaning and disinfection instructions. Any deviation from the manufacturer's instructions or agents listed below may impact the performance or durability of the product. See cleaning and disinfection manufacturer's instructions for additional warnings and cautions.
- See the **Instructions for Use** included with the product for additional warnings and cautions.

Warning

- Inspect the mask components after disinfection and replace any components that are worn or damaged. See the **Instructions for Use** included with the product for additional warnings and cautions.

Notes

- The headgear cannot be processed using the cleaning and disinfection methods listed. The headgear must be replaced or hand/machine washed and allowed to air dry before reuse, as indicated in the **Instructions for Use** included with the product.
- After disinfection, discoloration of the cushion and a slight odor is normal.

Cleaning

- Disassemble the mask per the **Instructions for Use** included with the product.
- Clean the mask using a soft bristle brush to adequately remove adhering substances from each component while soaking in a commercially available, anionic detergent such as MediZyme or Manu Klenz (used with chemical and gas plasma treatments). Extra attention should be given to the crevices and cavities during cleaning.

 NOTE: If compatible with the manufacturer's instructions, the detergent may be applied within the washer/disinfector cycle.

- Rinse the mask with potable water and air dry out of direct sunlight. If cleaning manually, use a minimum 5 liters per mask.

High-Level Disinfection

Mask and cushions may be reused after cleaning and disinfection by one of the chemical or gas plasma options below.

Chemical Disinfection: Apply only one treatment option following the manufacturer's instructions.

| TREATMENT | EXPOSURE TIME | ACTIVE INGREDIENTS | SOLUTION | RINSE | MAXIMUM CYCLES |
|-------------------------------------|---------------|--|----------------|--|----------------|
| Cidex Activated Dialdehyde Solution | 45 minutes | 2.4% Glutaraldehyde | Full Strength | Potable Water, 5 Liters per Mask (minimum) | 30 |
| or | | | | | |
| Cidex OPA | 12 minutes | 0.55% Ortho-Phthalaldehyde | Full Strength | Potable Water, 5 Liters per Mask (minimum) | 30 |
| or | | | | | |
| Wavicide | 45 minutes | 2.65% Glutaraldehyde | Full Strength | Potable Water, 5 Liters per Mask (minimum) | 30 |
| or | | | | | |
| Control III | 10 minutes | 10% Alkyl Dimethyl Benzyl Ammonium Chloride 10% Alkyl Dimethyl Ethyl Benzyl Ammonium Chloride | 8 ml per liter | Potable Water, 5 Liters per Mask (minimum) | 30 |

Continued

Hydrogen Peroxide Gas Plasma

| TREATMENT | EXPOSURE TIME | MAXIMUM CYCLES |
|--------------|--------------------|----------------|
| Sterrad 100S | Machine Controlled | 30 |

Visit us on-line at www.respironics.com or contact customer service to obtain mask cleaning and disinfection instructions.

Respironics: 1-800-345-6443 (USA or Canada only) or 724-387-4000

Respironics Deutschland: +49 (0) 8152 93060

Respironics Inc.
1001 Murry Ridge Lane
Murrysville, PA 15668
USA

Respironics Deutschland
Gewerbestr. 17
82211 Herrsching, Germany

